



CAREVENT POST USE TEST

Team	Date:
Team Member Number	

Demerit: 2 demerits for each omission

For Failure to:	Demerit
1 Remove and wash mask, connector and diaphragm	
2 Wipe resuscitator and hose with germicide or clean cloth	
3 Install full cylinders	
4 Lubricate o-ring	
5 Inspect exhalation valve	
6 Re-assemble resuscitator	
7 Connect facemask	
8 Perform field test	
9 Repack	

Total Demerits:

Signature of Judge: _____



ONTARIO MINE RESCUE
PREPARED SINCE 1929

CAREVENT POST USE TEST