



ONTARIO MINE RESCUE
PREPARED SINCE 1929

CARE-VENT POST USE SERVICE

<i>Team</i>	<i>Date:</i>
-------------	--------------

<i>Team Member:</i>

Demerit: 2 demerits for each omission, 1 demerit for wrong order

<i>For Failure to:</i>	<i>Demerit</i>
------------------------	----------------

1	Remove and wash mask, connector, and diaphragm	
---	--	--

2	Wipe resuscitator and hose with germicide or clean cloth	
---	--	--

3	Install full cylinders	
---	------------------------	--

4	Lubricate "o"-ring	
---	--------------------	--

5	Inspect exhalation valve	
---	--------------------------	--

6	Re-assemble resuscitator	
---	--------------------------	--

7	Connect facemask	
---	------------------	--

8	Perform field test	
---	--------------------	--

9	Repack	
---	--------	--

Total Demerits:

Signature of Judge
